




AJL BLUE



The future of vision

# AJL BLUE Ophthalmic Solution

|                                    |   |
|------------------------------------|---|
| <i>Brand name</i>                  | <b>AJL BLUE</b>   |
| <i>Reference</i>                   | <b>AJL BLUE 0.06%</b>   |
| <i>Description and indications</i> | <p>AJL BLUE 0.06% is an ophthalmic sterile trypan blue solution, intraocular and biocompatible used in ophthalmic surgery of the anterior segment. This ophthalmic dye is used as a coadjuvant during capsulorhexis and minimises the risk of laceration in eyes with mature cataracts, with poor reflexes at the back of the eye or with small pupils.</p> <ul style="list-style-type: none"> <li>• Capsulorhexis.</li> <li>• Phacoemulsification</li> <li>• Small incision cataract surgery.</li> <li>• Extracapsular cataract surgery.</li> <li>• Anterior/posterior capsulotomy.</li> <li>• Membrane peeling.</li> <li>• Improved visualization.</li> </ul> |
| <i>Composition</i>                 | Trypan blue sterile ophthalmic solution at 0.06%  |
| <i>Sterilization method</i>        | Steam   |
| <i>Supply</i>                      | 1 ml prefilled glass syringe and Luer lock 27G cannula  |
| <i>Shelf-life</i>                  | 3 years   |

\* Product with CE marking